



Typical Comparison Calibration


Additional Model SB148
10 channel Switchbox
(up to 60 channels)



Model SB148
10 channel Switchbox



This diagram shows a typical calibration system.
Call our application engineers to assist you in
the selection of equipment for your calibration
laboratory.




Model F700 Thermometry Bridge




Model RR
25 ohm resistor

Model RW
100 ohm resistor


FA-3
Junction Box



Ice Bath



Liquid Bath



Triple point
of water cell



PC with U-CAL
data logging software

IEEE

AUTOMATIC SYSTEMS LABORATORIES

275 King Henry's Drive, New Addington,
Croydon, CR0 0AE, UK.

Tel: +44 (0)1689 800799

Fax: +44(0)1689 800405

email: sales@aslltd.co.uk

web: www.aslltd.co.uk